



SHELLEY  
MEDICAL  
IMAGING  
TECHNOLOGIES

*Sales Office*

London, Ontario, Canada  
Phone: 1 (519) 690-0874  
Fax: 1 (519) 690-0875  
Email: [sales@simutech.com](mailto:sales@simutech.com)  
Web: [www.simutech.com](http://www.simutech.com)

WORLD LEADER  
IN ANATOMICAL HUMAN  
VASCULAR REPLICAS




## ABDOMINAL SOFT MEDIUM-LARGE ANEURYSM, LONG-LARGE NECK MODEL A-S-A-020

A smaller neck version is available; Model: A-S-A-006-B (shorter neck).



### Internal dimensions of Model A-S-A-020 (lumen):



These models are compatible with modern imaging modalities such as digital subtraction angiography, computed tomography and magnetic resonance imaging. Providing the use of an adequate circulating fluid, Doppler techniques can also be performed. The in vitro models transparency to light makes them suitable for video and photographic monitoring.