



SHELLEY  
MEDICAL  
IMAGING  
TECHNOLOGIES

*Sales Office*

London, Ontario, Canada

Phone: 1 (519) 690-0874

Fax: 1 (519) 690-0875


Email: sales@simutec.com

Web: www.simutec.com


**WORLD LEADER  
IN ANATOMICAL HUMAN  
VASUCULAR REPLICAS**




**HEAD H+N-R-A-002 MODEL**




Rigid Head model with its tubes



Tubes



Silicone block



These models are compatible with modern imaging modalities such as digital subtraction angiography, computed tomography and magnetic resonance imaging. Providing the use of an adequate circulating fluid, Doppler techniques can also be performed. The in vitro models transparency to light makes them suitable for video and photographic monitoring.