



SHELLEY
MEDICAL
IMAGING
TECHNOLOGIES

Sales Office

London, Ontario, Canada

Phone: 1 (519) 690-0874

Fax: 1 (519) 690-0875


Email: sales@simutec.com

Web: www.simutec.com


WORLD LEADER
IN ANATOMICAL HUMAN
VASUCULAR REPLICAS




THORACIC T-R-N-002 MODEL



Rigid thoracic model with its tubes



Tubes



Rigid silicone block



These models are compatible with modern imaging modalities such as digital subtraction angiography, computed tomography and magnetic resonance imaging. Providing the use of an adequate circulating fluid, Doppler techniques can also be performed. The in vitro models transparency to light makes them suitable for video and photographic monitoring.