



S H E L L E Y
M E D I C A L
I M A G I N G
T E C H N O L O G I E S

Sales Office

London, Ontario, Canada

Phone: 1 (519) 690-0874

Fax : 1 (519) 690-0875

Email: sales@simutec.com


Web: www.simutec.com

**WORLD LEADER
IN ANATOMICAL HUMAN
VASUCULAR REPLICAS**



HEART MODEL with RED AND BLUE DYE Model: T-S-N-020+ Complete

Left and right side filled with red and blue dye.



These models are compatible with modern imaging modalities such as digital subtraction angiography, computed tomography and magnetic resonance imaging. Providing the use of an adequate circulating fluid, Doppler techniques can also be performed. The in vitro models transparency to light makes them suitable for video and photographic monitoring.